



PeelVue+™ Sterilization Pouches and Sterility Shelf Life

Items packaged and sterilized in PeelVue+™ Sterilizing Pouches have an event-related shelf life. This means that the items within the PeelVue+™ packaging are considered sterile and remain sterile until an event occurs that breaks the sterile barrier. Events include tears, punctures, moisture, or opening of the package.

Because of the great variation in methods used to store sterile goods, it is not possible to make a general recommendation concerning the shelf life. However, based on logistics, maximum storage periods of 6 months have proven feasible provided that protected storage is ensured. This corresponds to the ISO 11607 recommendation (DIN 58 953 part 9 for the shelf-life of disposable items) as well as to the German RKI recommendation.

Sincerely,

Dux Dental B.V.

A handwritten signature in blue ink, consisting of a stylized 'J' followed by a horizontal line and a vertical stroke.

J. Doorman
Managing Director