

Safe, convenient, cost-effective and easy to use

Safety

External and **patented Internal Processing Indicators** built directly into each pouch.


Vents located on either side of the internal indicator to allow the same ambient conditions inside and outside the pouch.

10 mm broad, safe seal.

Convenience

A pre-applied sealing tape with a unique non-adhesive edge for easy peel and seal of the pouch – even with gloves.

The self-adhesive 19 mm wide sealing tape is designed into the pouch, making a heat sealing machine superfluous.



*Don't take safety for granted.
Use the PeelVue⁺ pouch.*



Various pouch sizes for all potential applications


PeelVue+ pouches are easy to peel and seal. They feature a patented Internal-Processing-Indicator (IPI), and are printed with lead-free inks on medical-grade, virgin, non-recycled paper. After removing the clear film, the paper side of the pouch can function as a sterilized tray liner.



Ref. number	Description		Size	Color Code	Quantity
31610	Implant, Dam Clamp		90 x 135 mm		200
31611	Handpiece, Hand Instrument	④	70 x 230 mm		200
31612	Handpiece, Hand Instrument	⑤	90 x 230 mm		200
31613	Burs	①	55 x 100 mm		200
31614	Impression Tray, Bur Block		135 x 165 mm		200
31615	Multiple Instruments and cassettes		215 x 380 mm		200
31616	Hand Instrument (long)		85 x 305 mm		200
31617	Surgical, Orthodontic		110 x 305 mm		200
31618	Multiple Instruments, and small cassettes	②	205 x 305 mm		200
31619	Multiple Instruments		135 x 305 mm		200
31643	Large cassettes	③	265 x 405 mm		200
31647	Large Cassettes (I.M.S. Type)		335 x 510 mm		200

The Smart Pouch...

...now with color identification!



PeelVue+ sterilization pouches comply with U.S. CDC, ADA and German RKI guidelines. ISO 11607